

CONDUCTIVITY METER WITH DIGITAL DISPLAY MODEL 103



SPECIFICATION

RANGE	0-20.00 mS/CM
RESOLUTION	10 μ S/CM
ACCURACY	$\pm 1\%$ ± 1 Digit
AUTOMATIC TEMPERATURE COMPENSATION	5-50 $^{\circ}$ C
CELL CONSTANT	1.0
CELL CONSTANT ADJUST	Yes
DISPLAY	12.7 mm (0.5") high LCD
POWER SOURCE	Internal 9 Volt battery
BATTERY LIFE	50 hours typical
DIMENSIONS (meter only)	137mm x 92mm x 48mm (5 3/8" x 3 5/8" x 1 7/8")
WEIGHT (meter and battery)	0.3 Kg (10.5 oz)

JENCO INSTRUMENTS, INC.

7968 ARJONS DRIVE, SUITE C,
SAN DIEGO, CA 92126
TLX: 324880 JENCOI
PHONE: (619) 578-2828
FAX: (619) 578-2886